

兔抗 FOXD1 多克隆抗体

中文名称：兔抗 FOXD1 多克隆抗体

英文名称：Anti-FOXD1 rabbit polyclonal antibody

别名：FKHL8, FREAC4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FOXD1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This intronless gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in tumor formation.
Applications:	ELISA, IHC
Name of antibody:	FOXD1
Immunogen:	Synthetic peptide of human FOXD1
Full name:	Forkhead box protein D1
Synonyms :	FKHL8, FREAC4
SwissProt:	Q16676

ELISA Recommended dilution:	500-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	10-50

