

兔抗 FOXD2 多克隆抗体

中文名称：兔抗 FOXD2 多克隆抗体

英文名称：Anti-FOXD2 rabbit polyclonal antibody

别名：FKHL17, FREAC9, FREAC-9

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FOXD2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. FOX (Forkhead box) proteins are a family of transcription factors that play important roles in regulating the expression of genes involved in cell growth, proliferation, differentiation, and longevity. Many FOX proteins are important to embryonic development.
Applications:	ELISA, IHC
Name of antibody:	FOXD2
Immunogen:	Synthetic peptide of human FOXD2
Full name:	forkhead box D2

Synonyms :	FKHL17, FREAC9, FREAC-9
SwissProt:	O60548
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	15-50

