

兔抗 FOXF1 多克隆抗体

中文名称：兔抗 FOXF1 多克隆抗体

英文名称：Anti-FOXF1 rabbit polyclonal antibody

别名：FKHL5; ACDMPV; FREAC1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FOXF1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in the regulation of pulmonary genes as well as embryonic development.
Applications:	ELISA, IHC
Name of antibody:	FOXF1
Immunogen:	Synthetic peptide of human FOXF1
Full name:	forkhead box F1
Synonyms :	FKHL5; ACDMPV; FREAC1

SwissProt:	Q12946
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human cervical cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

