

LMCD1 抗原（重组蛋白）

中文名称：LMCD1 抗原（重组蛋白）

英文名称：LMCD1 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 166-365 amino acids of human LMCD1

技术规格：

Full name:	LIM and cysteine rich domains 1
Swissprot:	Q9NZU5
Gene Accession:	BC000646
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the LIM-domain family of zinc finger proteins. The encoded protein contains an N-terminal cysteine-rich domain and two C-terminal LIM domains. The presence of LIM domains suggests involvement in protein-protein interactions. The protein may act as a co-regulator of transcription along with other transcription factors. Alternate splicing results in multiple transcript variants of this gene.