

## METTL1 抗原（重组蛋白）

中文名称：METTL1 抗原（重组蛋白）

英文名称：METTL1 Antigen (Recombinant Protein)

别名：TRM8; TRMT8; C12orf1; YDL201w

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-276 amino acids of human METTL1

技术规格：

<b>Full name:</b>	methyltransferase like 1
<b>Synonyms:</b>	TRM8; TRMT8; C12orf1; YDL201w
<b>Swissprot:</b>	Q9UBP6
<b>Gene Accession:</b>	BC000550
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene is similar in sequence to the <i>S. cerevisiae</i> YDL201w gene. The gene product contains a conserved S-adenosylmethionine-binding motif and is inactivated by phosphorylation. Alternative splice variants encoding different protein isoforms have been described for this gene. A pseudogene has been identified on chromosome X.