

兔抗 FOXO1(Phospho-Ser329) 多克隆抗体

中文名称：兔抗 FOXO1(Phospho-Ser329) 多克隆抗体

英文名称： Anti-FOXO1(Phospho-Ser329) rabbit polyclonal antibody

别名： FKH1; FKHR; FOXO1A

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： FOXO1(Phospho-Ser329)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in myogenic growth and differentiation. Translocation of this gene with PAX3 has been associated with alveolar rhabdomyosarcoma.
Applications:	WB
Name of antibody:	FOXO1(Phospho-Ser329)
Immunogen:	Peptide sequence around phosphorylation site of Serine 329 (R-L-S(p)-P-I) derived from Human FOX

	O1A.
Full name:	forkhead box O1
Synonyms :	FKH1; FKHR; FOXO1A
SwissProt:	Q12778
WB Predicted band size:	70 kDa
WB Positive control:	Hela cells lysates untreated or treated with serum
WB Recommended dilution:	500-1000

