

MRPL16 抗原（重组蛋白）

中文名称：MRPL16 抗原（重组蛋白）

英文名称：MRPL16 Antigen (Recombinant Protein)

别名：L16mt; MRP-L16; PNAS-111

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 30-251 amino acids of human MRPL16

技术规格：

Full name:	mitochondrial ribosomal protein L16
Synonyms:	L16mt; MRP-L16; PNAS-111
Swissprot:	Q9NX20
Gene Accession:	BC001040
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter c

ontain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein.