

MTHFD2L 抗原（重组蛋白）

中文名称：MTHFD2L 抗原（重组蛋白）

英文名称： MTHFD2L Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 136-268 amino acids of human MTHFD2L

技术规格：

Full name:	methylenetetrahydrofolate dehydrogenase (NADP+ dependent) 2 like
Swissprot:	Q9H903
Gene Accession:	BC032771
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The specific function of this protein remains unknown.