

NEGR1 抗原（重组蛋白）

中文名称：NEGR1 抗原（重组蛋白）

英文名称：NEGR1 Antigen (Recombinant Protein)

别名：Ntra; KILON; IGLON4; DMML2433

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 129-324 amino acids of human NEGR1

技术规格：

Full name:	neuronal growth regulator 1
Synonyms:	Ntra; KILON; IGLON4; DMML2433
Swissprot:	Q7Z3B1
Gene Accession:	BC036771
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in cell-adhesion. May function as a trans-neural growth-promoting factor in regenerative axon sprouting in the mammalian brain (By similarity).