

## 兔抗 FOXO4 多克隆抗体

- 中文名称：兔抗 FOXO4 多克隆抗体
- 英文名称：Anti-FOXO4 rabbit polyclonal antibody
- 别名：AFX; AFX1; MLLT7
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：FOXO4
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：Unconjugate



### 技术规格

#### Background:

This gene encodes a member of the O class of winged helix/forkhead transcription factor family. Proteins encoded by this class are regulated by factors involved in growth and differentiation indicating they play a role in these processes. A translocation involving this gene on chromosome X and the homolog of the Drosophila trithorax gene, encoding a DNA binding protein, located on chromosome 11 is associated with leukemia. Multiple transcript variants encoding different isoforms have been found for this gene.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FOXO4
<b>Immunogen:</b>	Synthetic peptide of human FOXO4
<b>Full name:</b>	forkhead box O4
<b>Synonyms :</b>	AFX; AFX1; MLLT7
<b>SwissProt:</b>	P98177
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and human tonsil
<b>IHC Recommend dilution:</b>	25-100





**mlbio** 梅联生物  
Good elisakit producers