

NID2 抗原（重组蛋白）

中文名称：NID2 抗原（重组蛋白）

英文名称：NID2 Antigen (Recombinant Protein)

别名：NID-2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1237-1375 amino acids of human NID2

技术规格：

Full name:	nidogen 2
Synonyms:	NID-2
Swissprot:	Q14112
Gene Accession:	BC035608
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the nidogen family of basement membrane proteins. This protein is a cell-adhesion protein that binds collagens I and IV and laminin and may be involved in maintaining the structure of the basement membrane.