

NIPBL 抗原(重组蛋白)

中文名称: NIPBL 抗原(重组蛋白)

英文名称: NIPBL Antigen (Recombinant Protein)

别 名: CDLS; IDN3; Scc2; CDLS1; IDN3-B

储 存: 冷冻(-20℃)

相关类别: 抗原

概述:

Fusion protein corresponding to a region derived from 2523-2687 amino acids of human NIPBL

技术规格:

Full name:	NIPBL, cohesin loading factor
Synonyms:	CDLS; IDN3; Scc2; CDLS1; IDN3-B
Swissprot:	Q6KC79
Gene Accession:	BC033847
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes the homolog of the Drosophila melanogast er Nipped-B gene product and fungal Scc2-type sister chromat id cohesion proteins. The Drosophila protein facilitates enhance r-promoter communication of remote enhancers and plays a rolle in developmental regulation. It is also homologous to a family of chromosomal adherins with broad roles in sister chromat



id cohesion, chromosome condensation, and DNA repair. The h uman protein has a bipartite nuclear targeting sequence and a putative HEAT repeat. Condensins, cohesins and other complex es with chromosome-related functions also contain HEAT repeats. Mutations in this gene result in Cornelia de Lange syndrom e, a disorder characterized by dysmorphic facial features, growth delay, limb reduction defects, and cognitive disability. Two transcript variants encoding different isoforms have been found for this gene.