

NKD2 抗原（重组蛋白）

中文名称：NKD2 抗原（重组蛋白）

英文名称：NKD2 Antigen (Recombinant Protein)

别名：Naked2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-311 amino acids of human NKD2

技术规格：

Full name:	naked cuticle homolog 2
Synonyms:	Naked2
Swissprot:	Q969F2
Gene Accession:	BC012176
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of a family of proteins that function as negative regulators of Wnt receptor signaling through interaction with Dishevelled family members. The encoded protein participates in the delivery of transforming growth factor alpha-containing vesicles to the cell membrane. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.