

## NUTF2 抗原（重组蛋白）

中文名称：NUTF2 抗原（重组蛋白）

英文名称：NUTF2 Antigen (Recombinant Protein)

别名：NTF2; PP15; NTF-2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-127 amino acids of human NUTF2

技术规格：

|                           |  |
|---------------------------|--|
| <b>Full name:</b>         | nuclear transport factor 2   |
| <b>Synonyms:</b>          | NTF2; PP15; NTF-2  |
| <b>Swissprot:</b>         | P61970   |
| <b>Gene Accession:</b>    | BC002348   |
| <b>Purity:</b>            | >85%, as determined by Coomassie blue stained SDS-PAGE   |
| <b>Expression system:</b> | Escherichia coli   |
| <b>Tags:</b>              | His tag C-Terminus, GST tag N-Terminus   |
| <b>Background:</b>        | This gene encodes a cytosolic factor that facilitates protein transport into the nucleus. The encoded protein is required for nuclear import of the small Ras-like GTPase, Ran which is involved in numerous cellular processes. This protein also interacts with the nuclear pore complex glycoprotein p62. |