

NXNL1 抗原(重组蛋白)

中文名称: NXNL1 抗原 (重组蛋白)

- 英文名称: NXNL1 Antigen (Recombinant Protein)
- 别名: nucleoredoxin like 1; RDCVF; TXNL6
- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述:

Fusion protein corresponding to a region derived from 1-212 amino acids of human NXNL1

技术规格:

Full name:	nucleoredoxin like 1
Synonyms:	RDCVF; TXNL6
Swissprot:	Q96CM4
Gene Accession:	BC014127
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Retinitis pigmentosa (RP) is a disease that leads to blindness by de generation of cone photoreceptors. Rods produce factors required f or cone viability. The protein encoded by this gene is one of those factors and is similar to a truncated form of thioredoxin. This gene has been proposed to have therapeutic value against RP.