

## 兔抗 FOXP4 多克隆抗体

中文名称：兔抗 FOXP4 多克隆抗体

英文名称： Anti-FOXP4 rabbit polyclonal antibody

别名： hFKHLA

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： FOXP4

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene belongs to subfamily P of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in the regulation of tissue- and cell type-specific gene transcription during both development and adulthood. Many members of the forkhead box gene family, including members of subfamily P, have roles in mammalian oncogenesis. This gene may play a role in the development of tumors of the kidney and larynx. Alternative splicing of this gene produces multiple transcript variants, some encoding different isoforms.

<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	FOXP4
<b>Immunogen:</b>	Synthetic peptide of human FOXP4
<b>Full name:</b>	forkhead box P4
<b>Synonyms :</b>	hFKHLA
<b>SwissProt:</b>	Q8IVH2
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	20-100
<b>WB Predicted band size:</b>	73 kDa
<b>WB Positive control:</b>	293T cell
<b>WB Recommended dilution:</b>	200-1000



