

PAF1 抗原（重组蛋白）

中文名称：PAF1 抗原（重组蛋白）

英文名称：PAF1 Antigen (Recombinant Protein)

别名：PAF1 homolog, Paf1/RNA polymerase II complex component; PD2; F23149_1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 332-531 amino acids of human PAF1

技术规格：

Full name:	PAF1 homolog, Paf1/RNA polymerase II complex component
Synonyms:	PD2; F23149_1
Swissprot:	Q8N7H5
Gene Accession:	BC000017
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a subunit of the polymerase associated factor (PAF1) complex. The PAF1 complex interacts with RNA polymerase II and plays a role in transcription elongation as well as histone modifications including ubiquitylation and methylation. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.