

兔抗 FRK 多克隆抗体

中文名称：兔抗 FRK 多克隆抗体

英文名称： Anti-FRK rabbit polyclonal antibody

别名： GTK; RAK; PTK5

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： FRK

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the TYR family of protein kinases. This tyrosine kinase is a nuclear protein and may function during G1 and S phase of the cell cycle and suppress growth.
Applications:	ELISA, WB, IHC
Name of antibody:	FRK
Immunogen:	Fusion protein of human FRK
Full name:	fyn-related kinase
Synonyms :	GTK; RAK; PTK5
SwissProt:	P42685
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human colon cancer and Human brain
IHC Recommend dilution:	50-200
WB Predicted band size:	58 kDa
WB Positive control:	A431 cells
WB Recommended dilution:	500-2000



