

兔抗 FSCB 多克隆抗体

中文名称: 兔抗 FSCB 多克隆抗体

英文名称: Anti-FSCB rabbit polyclonal antibody

别名: fibrous sheath CABYR binding protein; C14orf155

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: FSCB

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | May be involved in the later stages of fibrous sheath biogenesis and spermatozoa capacitation. Inhibits ROPN1 and ROPN1L SUMOylation. Binds calcium. |
| Applications: | ELISA, IHC |
| Name of antibody: | FSCB |
| Immunogen: | Fusion protein of human FSCB |
| Full name: | fibrous sheath CABYR binding protein |
| Synonyms: | C14orf155 |
| SwissProt: | Q5H9T9 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human esophagus cancer and Human colorectal ca |

| | |
|--------------------------------|--------|
| | ncer |
| IHC Recommend dilution: | 50-200 |



