

PER3 抗原(重组蛋白)

中文名称: PER3 抗原(重组蛋白)

英文名称: PER3 Antigen (Recombinant Protein)

别名: GIG13; FASPS3

储 存: 冷冻(-20℃)

相关类别: 抗原

概述:

Fusion protein corresponding to a region derived from 179-378 amino acids of human PER3

技术规格:

Full name:	period circadian regulator 3
Synonyms:	GIG13; FASPS3
Swissprot:	P56645
Gene Accession:	BC026102
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is a member of the Period family of genes and is expres sed in a circadian pattern in the suprachiasmatic nucleus, the prima ry circadian pacemaker in the mammalian brain. Genes in this famil y encode components of the circadian rhythms of locomotor activit y, metabolism, and behavior. This gene is upregulated by CLOCK/A RNTL heterodimers but then represses this upregulation in a feedba ck loop using PER/CRY heterodimers to interact with CLOCK/ARNTL.



Polymorphisms in this gene have been linked to sleep disorders. M
ultiple transcript variants encoding different isoforms have been fou
nd for this gene.