

## PFDN2 抗原（重组蛋白）

中文名称：PFDN2 抗原（重组蛋白）

英文名称：PFDN2 Antigen (Recombinant Protein)

别名：PFD2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-154 amino acids of human PFDN2

技术规格：

<b>Full name:</b>	prefoldin subunit 2
<b>Synonyms:</b>	PFD2
<b>Swissprot:</b>	Q9UHV9
<b>Gene Accession:</b>	BC012464
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the prefoldin beta subunit family. The encoded protein is one of six subunits of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils.