

## 小鼠抗 mRaspberry 单克隆抗体

中文名称: 小鼠抗 mRaspberry 单克隆抗体

英文名称: Anti-mRaspberry mouse monoclonal antibody

相关类别: 标签抗体

储 存: 冷冻 (-20℃) 避光

宿 主: Mouse

抗 原: mRaspberry

反应种属: mRaspberry, AmCyan1

标 记 物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	NULL
<b>Applications:</b>	WB
<b>Name of antibody:</b>	mRaspberry
<b>Immunogen:</b>	Purified mRaspberry protein expressed in HEK293 cell
<b>WB Predicted band size:</b>	26 kDa
<b>WB Positive control:</b>	HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY mRaspberry ne) cDNA for 48 hrs and lysed
<b>WB Recommended dilution:</b>	1000-5000