

## 小鼠抗 ZsGreen1 单克隆抗体

中文名称： 小鼠抗 ZsGreen1 单克隆抗体

英文名称： Anti-ZsGreen1 mouse monoclonal antibody

相关类别： 标签抗体

储 存： 冷冻（-20℃） 避光

宿 主： Mouse

抗 原： ZsGreen1

反应种属： ZsGreen1

标 记 物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	ZsGreen1 is a bright green fluorescent protein derived from a Zoanthus sp. reef coral that has been modified for high solubility, bright emission, and rapid chromophore maturation. The ZsGreen1 fluorescent protein has been reported to have better expression than EGFP in lung and muscle tissue and in cultured cells.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ZsGreen1
<b>Immunogen:</b>	Purified ZsGreen1 protein expressed in HEK293 cell
<b>WB Predicted band size:</b>	26 kDa
<b>WB Positive control:</b>	HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZsGreen1

	ne) cDNA for 48 hrs and lysed
<b>WB Recommended dilution:</b>	1000-5000