

小鼠抗 ACTB 单克隆抗体

中文名称： 小鼠抗 ACTB 单克隆抗体

英文名称： Anti-ACTB mouse monoclonal antibody

别名： Beta actin

相关类别： 内参抗体

储存： 冷冻（-20℃）

宿主： Mouse

抗原： ACTB

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.
Applications:	WB
Name of antibody:	ACTB
Immunogen:	Fusion protein of human ACTB
Full name:	Beta actin
SwissProt:	P60709
WB Predicted band size:	42 kDa

WB Positive control:	Hela cell, 293T cell, Human cerebella tissue, Mouse brain tissue, Rat brain tissue, NIH/3T3 cell, Jurkat cell lysates
WB Recommended dilution:	1000-5000

