

兔抗 TMSB4X 多克隆抗体

- 中文名称：兔抗 TMSB4X 多克隆抗体
- 英文名称：Anti-TMSB4X rabbit polyclonal antibody
- 别名：FX; TB4X; PTMB4; TMSB4
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：TMSB4X
- 反应种属：Human, Mouse, Rat
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an actin sequestering protein which plays a role in regulation of actin polymerization. The protein is also involved in cell proliferation, migration, and differentiation. This gene escapes X inactivation and has a homolog on chromosome Y.
Applications:	ELISA, IHC
Name of antibody:	TMSB4X
Immunogen:	Synthetic peptide of human TMSB4X
Full name:	thymosin beta 4, X-linked
Synonyms :	FX; TB4X; PTMB4; TMSB4
SwissProt:	P62328

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100

