

兔抗 TMSB10 多克隆抗体

中文名称：兔抗 TMSB10 多克隆抗体

英文名称：Anti-TMSB10 rabbit polyclonal antibody

别名：TB10; MIG12

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TMSB10

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Thymosin beta-10 plays an important role in the organization of the cytoskeleton. Thymosin beta-10 binds to and sequesters actin monomers (G actin) and therefore inhibits actin polymerization. The expression of thymosin beta-10 dramatically decreases after birth.
Applications:	ELISA, IHC
Name of antibody:	TMSB10
Immunogen:	Synthetic peptide of human TMSB10
Full name:	thymosin beta 10
Synonyms :	TB10; MIG12
SwissProt:	P63313

ELISA Recommended dilution:	2000-10000
IHC positive control:	Human breast cancer and Human brain
IHC Recommend dilution:	30-150

