

小鼠抗 METTL16 单克隆抗体

- 中文名称: 小鼠抗 METTL16 单克隆抗体
- 英文名称: Anti-METTL16 mouse monoclonal antibody
- 别 名: METT10D
- 储 存:冷冻(-20℃) 避光
- 抗 原: METTL16
- 宿 主: Mouse
- 反应种属: Human, Monkey
- 相关类别: 一抗
- 标记物: Unconjugate
- 克隆类型: mouse monoclonal

技术规格

	METT10D (methyltransferase 10 domain containing), a
	lso known as METTL16 (methyltransferase-like protein
	16) is a 562 amino acid protein that belongs to the
	METTL16/rlmF family and methyltransferase superfamil
	y. A putative methyltransferase, METT10D exists as tw
Background:	o alternatively spliced isoforms and undergoes phosp
	horylation in response to DNA damage, most likely b
	y Atm or ATR. The gene encoding METT10D maps to
	human chromosome 17, which comprises over 2.5% o
	f the human genome and encodes over 1,200 genes.
	Two key tumor suppressor genes are associated with



	chromosome 17, namely, p53 and BRCA1. Malfunctio n or loss of p53 expression is associated with malign ant cell growth and Li-Fraumeni syndrome. Like p53, BRCA1 is directly involved in DNA repair, and is linke d to predisposition of cancers of the ovary, colon, pr ostate gland and fallopian tubes.
Applications:	WB, IHC, IF
Name of antibody:	METTL16
Immunogen:	Fusion protein of human METTL16
Full name:	methyltransferase like 16 (METTL16)
Synonyms:	METT10D
SwissProt:	Q86W50
IHC positive control:	adenocarcinoma of human colon tissue and human p ancreas tissue; carcinoma of human prostate tissue a nd carcinoma of human bladder tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	64 kDa
WB Positive control:	HepG2, Hela, A549, COS7, Jurkat, MCF-7 cell lysates
WB Recommended dilution:	500-2000