

## 小鼠抗 MLF1 单克隆抗体

中文名称: 小鼠抗 MLF1 单克隆抗体

英文名称: Anti-MLF1 mouse monoclonal antibody

别名: myeloid leukemia factor 1 variant 1; myeloid leukemia factor 1 variant 2; myeloid leukemia factor 1 variant 3; myeloid leukemia fact

储存: 冷冻 (-20℃) 避光

抗原: MLF1

宿主: Mouse

反应种属: Human , Monkey

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoietic cells. Translocations between this gene and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq].
<b>Applications:</b>	WB, IF, FC
<b>Name of antibody:</b>	MLF1
<b>Immunogen:</b>	Fusion protein of human MLF1

<b>Full name:</b>	myeloid leukemia factor 1 (MLF1), transcript variant 1
<b>Synonyms:</b>	myeloid leukemia factor 1 variant 1; myeloid leukemia factor 1 variant 2; myeloid leukemia factor 1 variant 3; myeloid leukemia factor 1
<b>SwissProt:</b>	P58340
<b>WB Predicted band size:</b>	31 kDa
<b>WB Positive control:</b>	COS7 cell lysate
<b>WB Recommended dilution:</b>	200-1000