

小鼠抗 MYL4 单克隆抗体

中文名称： 小鼠抗 MYL4 单克隆抗体

英文名称： Anti-MYL4 mouse monoclonal antibody

别名： ALC1; AMLC; GT1; PRO1957

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： MYL4

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	Myosin is a hexameric ATPase cellular motor protein. It is composed of two myosin heavy chains, two nonphosphorylatable myosin alkali light chains, and two phosphorylatable myosin regulatory light chains. This gene encodes a myosin alkali light chain that is found in embryonic muscle and adult atria. Two alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008].
Applications:	WB, IHC
Name of antibody:	MYL4
Immunogen:	Fusion protein of human MYL4

Full name:	myosin, light chain 4, alkali; atrial, embryonic (MYL4), transcript variant 2
Synonyms:	ALC1; AMLC; GT1; PRO1957
SwissProt:	P12829
IHC positive control:	human adult heart tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	22 kDa
WB Positive control:	K562 cell lysate
WB Recommended dilution:	500-2000