

## 小鼠抗 MYL6B 单克隆抗体

中文名称：小鼠抗 MYL6B 单克隆抗体

英文名称：Anti-MYL6B mouse monoclonal antibody

别名：myosin light chain 6B; MLC1SA

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：MYL6B

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：mouse monoclonal

### 技术规格

<b>Background:</b>	Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain expressed in both slow-twitch skeletal muscle and in nonmuscle tissue. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MYL6B
<b>Immunogen:</b>	Fusion protein of human MYL6B
<b>Full name:</b>	myosin light chain 6B

<b>Synonyms:</b>	MLC1SA
<b>SwissProt:</b>	P14649
<b>WB Predicted band size:</b>	23 KD
<b>WB Positive control:</b>	K562 cell, Mouse testis tissue lysates
<b>WB Recommended dilution:</b>	200-1000