

小鼠抗 NFKBIA 单克隆抗体

中文名称: 小鼠抗 NFKBIA 单克隆抗体

英文名称: Anti-NFKBIA mouse monoclonal antibody

别 名: IKBA; MAD-3; NFKBI

储 存: 冷冻(-20℃) 避光

抗 原: NFKBIA

宿 主: Mouse

反应种属: Human

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a member of the NF-kappa-B inhibit or family, which contain multiple ankrin repeat domains. The encoded protein interacts with REL dimers to inhibit NF-kappa-B/REL complexes which are involved in inflam matory responses. The encoded protein moves between the cytoplasm and the nucleus via a nuclear localization signal and CRM1-mediated nuclear export. Mutations in this gene have been found in ectodermal dysplasia anhi drotic with T-cell immunodeficiency autosomal dominant disease. [provided by RefSeq, Aug 2011].
Applications:	WB, IF



Name of antibody:	NFKBIA
Immunogen:	Fusion protein of human NFKBIA
Full name:	nuclear factor of kappa light polypeptide gene enhancer
	in B-cells inhibitor, alpha (NFKBIA)
Synonyms:	IKBA; MAD-3; NFKBI
SwissProt:	P25963
WB Predicted band size:	36 kDa
WB Positive control:	HepG2, Hela, SVT2, A549, COS7 cell lysates
WB Recommended dilution:	500-2000