

小鼠抗 SAFB 单克隆抗体

中文名称: 小鼠抗 SAFB 单克隆抗体

英文名称: Anti-SAFB mouse monoclonal antibody

别 名: scaffold attachment factor B; HAP; HET; SAF-B; SAFB1; SAB-B1; SAF-B1

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Mouse

抗 原: SAFB

反应种属: Human

标 记 物: Unconjugate

克隆类型: Mouse Monoclonal

技术规格

Background:

This gene encodes a DNA-binding protein which has high s pecificity for scaffold or matrix attachment region DNA elem ents (S/MAR DNA). This protein is thought to be involved in attaching the base of chromatin loops to the nuclear matrix but there is conflicting evidence as to whether this protein is a component of chromatin or a nuclear matrix protein. Scaffold attachment factors are a specific subset of nuclear matrix proteins (NMP) that specifically bind to S/MAR. The encoded protein is thought to serve as a molecular base to assemble a 'transcriptosome complex' in the vicinity of actively transcribed genes. It is involved in the regulation of hea



	t shock protein 27 transcription, can act as an estrogen rece ptor co-repressor and is a candidate for breast tumorigenesi s. This gene is arranged head-to-head with a similar gene w hose product has the same functions. Multiple transcript var iants encoding different isoforms have been found for this gene.
Applications:	WB, IF
Name of antibody:	SAFB
Immunogen:	Fusion protein of human SAFB
Full name:	scaffold attachment factor B
Synonyms:	HAP; HET; SAF-B; SAFB1; SAB-B1; SAF-B1
SwissProt:	Q15424
WB Predicted band size:	103 kDa
WB Positive control:	A431, MDA-MB-468, Breast cancer, ZR75-1, MCF7 and Hela NE cell lysates
WB Recommended dilution:	1000-5000
IF positive control:	Hela cells
IF Recommend dilution:	100-200





