

兔抗 TOP1MT 多克隆抗体

中文名称：兔抗 TOP1MT 多克隆抗体

英文名称： Anti-TOP1MT rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻 (-20℃)

宿 主： Rabbit

抗 原： TOP1MT

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a mitochondrial DNA topoisomerase that plays a role in the modification of DNA topology. The encoded protein is a type IB topoisomerase and catalyzes the transient breaking and rejoining of DNA to relieve tension and DNA supercoiling generated in the mitochondrial genome during replication and transcription. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, IHC
Name of antibody:	TOP1MT
Immunogen:	Synthetic peptide of human TOP1MT
Full name:	topoisomerase (DNA) I, mitochondrial
SwissProt:	Q969P6

IHC positive control:	Human gastric cancer and human lung cancer
IHC Recommend dilution:	42510

