

小鼠抗 SMYD4 单克隆抗体

中文名称： 小鼠抗 SMYD4 单克隆抗体

英文名称： Anti-SMYD4 mouse monoclonal antibody

别名： SET and MYND domain containing 4; ZMYND21

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： SMYD4

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:

The SMYD (SET and MYND domain-containing) family of proteins, designated SMYD1, SMYD2, SMYD3 and SMYD4, play a role in transcriptional regulation and may interact with the HDAC family of histone deacetylases. Members of the SMYD family contain MYND (Myeloid translocation protein 8 (MTG8/ETO), Nerve protein and Deaf-1) and SET domains through which protein-protein interactions and transcriptional-associated activities are conferred. SMYD4 (SET and MYND domain-containing protein 4), also known as ZMYND21, is an 804 amino acid protein that contains one SET domain and one MYND-type zinc finger.

Applications:	WB
Name of antibody:	SMYD4
Immunogen:	Fusion protein of human SMYD4
Full name:	SET and MYND domain containing 4
Synonyms:	ZMYND21
SwissProt:	Q8IYR2
WB Predicted band size:	89 kDa
WB Positive control:	A431, MCF7 and Hela cell lysates
WB Recommended dilution:	200-1000

