

小鼠抗 SOCS4 单克隆抗体

中文名称： 小鼠抗 SOCS4 单克隆抗体

英文名称： Anti-SOCS4 mouse monoclonal antibody

别名： suppressor of cytokine signaling 4; SOCS7

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： SOCS4

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene contains a SH2 domain and a SOCS BOX domain. The protein thus belongs to the suppressor of cytokine signaling (SOCS), also known as STAT-induced STAT inhibitor (SSI), protein family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling. Two alternatively spliced transcript variants encoding the same protein have been found for this gene.
Applications:	WB
Name of antibody:	SOCS4
Immunogen:	Fusion protein of human SOCS4

Full name:	suppressor of cytokine signaling 4
Synonyms:	SOCS7
SwissProt:	Q8WXH5
WB Predicted band size:	51 KD
WB Positive control:	Hela cell lysate
WB Recommended dilution:	200-1000