

小鼠抗 TFR2 单克隆抗体

中文名称: 小鼠抗 TFR2 单克隆抗体

英文名称: Anti-TFR2 mouse monoclonal antibody

别 名: HFE3; TFRC2

抗 原: TFR2

储 存:冷冻(-20℃) 避光

宿 主: Mouse

相关类别: 一抗

反应种属: Human, Rat, Monkey, Mouse

标 记 物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a single-pass type II membrane protein, which is a member of the transferrin receptor-like family. This protein m ediates cellular uptake of transferrin-bound iron, and may be involved in iron metabolism, hepatocyte function and erythrocyte differ entiation. Mutations in this gene have been associated with hereditary hemochromatosis type III. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq].
Applications:	WB, IF, FC
Name of antibody:	TFR2
Immunogen:	Fusion protein of human TFR2



Full name:	transferrin receptor 2 (TFR2), transcript variant 1
Synonyms:	HFE3; TFRC2
SwissProt:	Q9UP52
WB Predicted band si ze:	89 kDa
WB Positive control:	HepG2, Hela, SVT2, A549, COS7, MDCK, PC12, MCF-7 cell lysates
WB Recommended dil ution:	500-2000