

## 小鼠抗 TF 单克隆抗体

- 中文名称: 小鼠抗 TF 单克隆抗体
- 英文名称: Anti-TF mouse monoclonal antibody
- 别名: transferrin; TFQTL1; PRO1557; PRO2086; HEL-S-71p
- 相关类别: 一抗
- 储存: 冷冻(-20℃)
- 宿 主: Mouse
- 抗 原: TF
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: mouse monoclonal

## 技术规格



Name of antibody:	TF
Immunogen:	Fusion protein of human TF
Full name:	transferrin
Synonyms:	TFQTL1; PRO1557; PRO2086; HEL-S-71p
SwissProt:	P02787
IHC positive control:	Human liver cancer
IHC Recommend diluti on:	100-500
WB Predicted band siz e:	77 kDa
WB Positive control:	Human tonsil tissue lysate
WB Recommended dil ution:	500-2000





