

小鼠抗 USP36 单克隆抗体

中文名称: 小鼠抗 USP36 单克隆抗体

英文名称: Anti-USP36 mouse monoclonal antibody

别名: DUB1

抗原: USP36

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Dog

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	Modification of cellular proteins by ubiquitin is an essential regulatory mechanism controlled by the coordinated action of multiple ubiquitin-conjugating and deubiquitinating enzymes. USP36 belongs to a large family of cysteine proteases that function as deubiquitinating enzymes (Quesada et al., 2004 [PubMed 14715245]).[supplied by OMIM].
Applications:	WB, IHC, IF
Name of antibody:	USP36
Immunogen:	Fusion protein of human USP36
Full name:	ubiquitin specific peptidase 36 (USP36)
Synonyms:	DUB1
SwissProt:	Q9P275
IHC positive control:	carcinoma of human lung tissue

IHC Recommend dilution:	30-150
WB Predicted band size:	123 kDa
WB Positive control:	MDCK cell lysate
WB Recommended dilution:	200-1000