

## SNX12 抗原（重组蛋白）

中文名称： SNX12 抗原（重组蛋白）

英文名称： SNX12 Antigen (Recombinant Protein)

别名： sorting nexin 12

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to N terminal 122 amino acids of human SNX12

### 技术规格

<b>Full name:</b>	sorting nexin 12
<b>Swissprot:</b>	Q9UMY4
<b>Gene Accession:</b>	BC103848
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the sorting nexin family. Members of this family contain a phosphoinositide binding domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members. A similar protein in mouse may be involved in regulating the neurite outgrowth. Alternate splicing results in multiple transcript variants.