

VEGFA 抗原

中文名称: VEGFA 抗原

英文名称: VEGFA Antigen

相关类别: 抗原

储 存: 冷冻(-20℃)

概述

Synthetic peptide corresponding to a region derived from 130-143 amino acids of human vascular endothelial growth factor A

技术规格

Full name	vascular endothelial growth factor A
Synonyms	VPF, VEGF, MVCD1
Swissprot	P15692
Gene Accession	NP_001020537
Background	This gene is a member of the PDGF/VEGF growth factor family and encodes a protein that is often found as a disulfide linked homodimer. This protein is a glycosylated mitogen that specific ally acts on endothelial cells and has various effects, including mediating increased vascular permeability, inducing angiogenesi s, vasculogenesis and endothelial cell growth, promoting cell m igration, and inhibiting apoptosis. Elevated levels of this protei n is linked to POEMS syndrome, also known as Crow-Fukase s yndrome. Mutations in this gene have been associated with pr oliferative and nonproliferative diabetic retinopathy. Alternativel y spliced transcript variants, encoding either freely secreted or



cell-associated isoforms, have been characterized. There is also evidence for the use of non-AUG (CUG) translation initiation si tes upstream of, and in-frame with the first AUG, leading to a dditional isoforms.