

ZNF512 抗原（重组蛋白）

中文名称： ZNF512 抗原（重组蛋白）

英文名称： ZNF512 Antigen (Recombinant Protein)

别名： zinc finger protein 512

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 368-567 amino acids of human ZNF512

技术规格

Full name:	zinc finger protein 512
Swissprot:	Q96ME7
Gene Accession:	BC043221
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein containing four putative zinc finger motifs. Zinc finger motifs may bind to proteins or nucleic acids. Zinc finger-containing proteins are involved in a variety of processes, including regulation of transcription. Alternative splicing results in multiple transcript variants for this gene.

mlbio 梅联生物
Good elisakit producers