

牛 IgG (H+L)

中文名称: 牛 IgG (H+L)

英文名称: Bovine IgG (H+L)

相关类别: IgG

储 存: 冷冻(-20℃) 避光

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant protein s in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defense against microorganis ms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions, to bind to the pathogen that elicited the response and to recruit other cells and molecules to dest roy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 v ariants.
Applications	Optimal tested applications and working dilutions must be d etermined by the end user.