

## 兔抗 CD14 多克隆抗体

中文名称：兔抗 CD14 多克隆抗体

英文名称：Anti-CD14 rabbit polyclonal antibody

相关类别：一抗

储存：冷冻（-20℃）避光

抗原：CD14

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a surface antigen that is preferentially expressed on monocytes/macrophages. It cooperates with other proteins to mediate the innate immune response to bacterial lipopolysaccharide. Alternative splicing results in multiple transcript variants encoding the same protein.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	CD14
<b>Immunogen:</b>	Fusion protein of human CD14
<b>Full name:</b>	CD14 molecule
<b>SwissProt:</b>	P08571
<b>IHC positive control:</b>	Human thyroid cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	40 kDa

WB Positive control:	A549 and NIH/3T3 cell
WB Recommended dilution:	500-2000



