

兔抗 CD28 多克隆抗体

中文名称：兔抗 CD28 多克隆抗体

英文名称：Anti-CD28 rabbit polyclonal antibody

别名：CD28 molecule; Tp44

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CD28

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|------------------------------------|---|
| Background: | The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. |
| Applications: | ELISA, IHC |
| Name of antibody: | CD28 |
| Immunogen: | Synthetic peptide of human CD28 |
| Full name: | CD28 molecule |
| Synonyms: | Tp44 |
| SwissProt: | P10747 |
| ELISA Recommended dilution: | 5000-10000 |

| | |
|--------------------------------|--------------|
| IHC positive control: | Human tonsil |
| IHC Recommend dilution: | 10-50 |

