

兔抗 DOCK1 多克隆抗体

中文名称: 兔抗 DOCK1 多克隆抗体

英文名称: Anti-DOCK1 rabbit polyclonal antibody

别名: ced5; DOCK180

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: DOCK1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene product binds to the SH3 domain of CRK protein. It may regulate cell surface extension and may have a role in the cell surface extension of an engulfing cell around a dying cell during apoptosis.
Applications:	ELISA, IHC
Name of antibody:	DOCK1
Immunogen:	Synthetic peptide of human DOCK1
Full name:	dedicator of cytokinesis 1
Synonyms :	ced5; DOCK180
SwissProt:	Q14185
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human gastric cancer

IHC Recommend dilution:

25-100



