

## 兔抗 DPH7 多克隆抗体

中文名称: 兔抗 DPH7 多克隆抗体

英文名称: Anti-DPH7 rabbit polyclonal antibody

别名: RRT2; WDR85; C9orf112

抗原: DPH7

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	WDR85 is a WD repeat-containing protein that plays a role in the first step of diphthamide biosynthesis (Carette et al., 2009 [PubMed 19965467]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	DPH7
<b>Immunogen:</b>	Fusion protein of human DPH7
<b>Full name:</b>	diphthamide biosynthesis 7
<b>Synonyms :</b>	RRT2; WDR85; C9orf112
<b>SwissProt:</b>	Q9BTV6
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

