

订购热线: 4008-898-798

Anti-NEGR1 antibody

Cat. No. ml125672

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-NEGR1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human NEGR1

Reactivity Human, Mouse, Rat

Content 0.78 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol NEGR1

Full name neuronal growth regulator 1

Synonyms Ntra; KILON; IGLON4; DMML2433

Swissprot Q7Z3B1

Target Background

May be involved in cell-adhesion. May function as a trans-neural growth-promoting factor in regenerative axon sprouting in the mammalian brain (By similarity).




订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human tonsil Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml125672(NEGR1 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000


联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn

Predicted cell location: Cytoplasm Positive control: Human gastric cancer

Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using ml125672(NEGR1 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)