

Anti-SULT4A1 antibody

 Cat. No.
 ml125693

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-SULT4A1 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human SULT4A1
Reactivity	Human, Mouse, Rat
Content	0.9 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	SULT4A1
Full name	sulfotransferase family 4A member 1

sulfotransferase family 4A member 1 NST; BRSTL1; SULTX3; BR-STL-1; DJ388M5.3; hBR-STL-1 Q9BR01

Target Background


Synonyms Swissprot

This gene encodes a member of the sulfotransferase family. The encoded protein is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. Polymorphisms in this gene may be associated with susceptibility to schizophrenia.



订购热线: 4008-898-798

Applications Immunohistochemistry Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml125693(SULT4A1 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798,021-61725725

联系QQ: 2881505695,2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn