

订购热线: 4008-898-798

Anti-RNF13 antibody

Cat. No. ml125825

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-RNF13 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human RNF13

Reactivity Human, Mouse, Rat

Content 1.08 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol RNF13

Full name ring finger protein 13

Synonyms RZF Swissprot O43567

Target Background

The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. The specific function of this gene has not yet been determined. Alternatively spliced transcript variants that encode the same protein have been reported. A pseudogene, which is also located on chromosome 3, has been defined for this gene.





订购热线: 4008-898-798

Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 50-300





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml125825(RNF13 Antibody) at dilution 1/60, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn