

Anti-SHC4 antibody

Cat. No.	ml125947
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-SHC4 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human SHC4
Reactivity	Human
Content	1.56 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	SHC4
Full name	SHC adaptor protein 4
Synonyms	RaLP; SHCD
Swissprot	Q6S5L8

Target Background

Activates both Ras-dependent and Ras-independent migratory pathways in melanomas. Contributes to the early phases of agrin-induced tyrosine phosphorylation of CHRNB1.

订购热线: 4008-898-798



Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human liver cancer

Recommended dilution: 90-450





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml125947(SHC4 Antibody) at dilution 1/90, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Cytoplasm

Positive control: Human cervical cancer

Recommended dilution: 90-450



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml125947(SHC4 Antibody) at dilution 1/90, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn